

# **EXERCISE PUTTY**

# Thera Putty

### **EXERCISES**

Below are a few exercises that are possible with putty. Using a higher level resistance putty, larger quantity of putty, and/or increasing number of repetitions will make exercise more challenging. Putty can also be used to exercise toes!

Levels of resistance from easiest to most difficult: Tan - Yellow - Red - Green - Blue - Black

#### **WARNING**

Therapy routines should be determined by a doctor, therapist, coach, or trainer.

### HAND FLEXION

**FINGER EXTENSION** 

putty held over the fingertip.

**FINGER SCISSOR** 

shape putty and repeat.

Place putty in the palm at the base of the fingers. Squeeze with fingers into the putty like you are making a fist. Release fingers, reshape putty and repeat process.

Roll putty into cylinder and wrap putty

around bent finger. Hold the putty with the other hand and try extending bent finger. The degree of resistance

is controlled by the thickness of the

Roll putty into a ball shape. Place be-


tween two fingers and squeeze them

together in a scissor-like motion. Re-



## THUMB PRESS

Roll putty into a ball and place in the palm of your hand. Press thumb into the putty towards small finger. Reshape putty and repeat.



### THUMB EXTENSION

Wrap putty around the tip of the thumb while it is bent. Hold the loose ends and try to straighten (extend) your thumb.



### THUMB PINCH

Make fist and put ball of putty between thumb and index finger. Pinch putty with thumb towards inside of



index finger.



### THUMB ADDUCTION

With hand opened, press the putty with your thumb against your index finger. Reshape putty and repeat exercise.




### **FINGER GRIP**

Make hook with fingertips and place putty ball between fingertips and palm. Press fingertips into putty.




### SCISSOR SPREAD

Wrap putty around tops of two fingers while they are together (steps 1 and 2). Spread fingers apart (step 3).



### **FINGER SPREAD**

Form putty into a pancake shape. Place on table or mold over finger tips while fingers are bunched together. Try to spread fingers apart.



### **Fabrication Enterprises Inc.**

PO Box 1500 White Plains, NY 10602 (USA) tel: +1-914-345-9300 • 800-431-2830 fax: +1-914-345-9800 • 800-634-5370 FabEnt.com







Authorized CE representative: RMS UK Ltd. 28 Trinity Road Nailsea, Somerset BS48 4NU (UK)

### **FINGER PINCH**

Roll putty into a ball shape. Pinch between the thumb and fingertips. This exercise can also be performed with the thumb and each finger separately. Reshape putty and repeat.



