

**Acti**move®



*Actimove® Manus Air*

Our first ready-to-wear stabilizing fracture wrist brace

# What If?

Distal radius fractures (DRFs) account for approximately one-sixth of all fractures. With current solutions, patients can experience discomfort, limited hand motion and difficulty to perform personal hygiene.<sup>1</sup> Healthcare professionals, in the meantime, strive to deliver the best possible outcomes for patients while battling the constant pressures of time and workload.

But what if there was a treatment that could maintain quality of care, while saving time? What if there was a solution that provides sufficient stabilization while still enabling patients to continue with activities of daily living?



## What if? Becomes Reality

Introducing Actimove® Manus Air. A patent-pending, ready-to-wear stabilizing fracture wrist brace that is up to 50% quicker to apply than a circumferential thermoplastic splint, water-resistant for improved hand hygiene and designed to allow easy post-op wound inspection.

<sup>1</sup>Source: Trends in the Incidence and Treatment of Distal Radius Fractures in the United States in Privately Insured and Medicare Advantage Enrollees Kakar S. et al 2022 Sanjeev Kakar<sup>1</sup>, Mohamed Nouredin<sup>2</sup>, Holly K. Van Houten<sup>1,3</sup>, Raphael Mwangi<sup>1</sup>, and Lindsey R. Sangaralingham<sup>1,3</sup>



[Full Study](#)



Detachable  
wrist strap  
for enhanced  
stability



# Reasons to adopt Actimove Manus Air

## Effective

- Supports fracture healing by immobilizing the wrist through the anatomically-shaped aluminium splint

## Time-Saving

- Faster alternative to thermoplastic splints – can be applied and removed within minutes

## Convenient

- Allows for easy wound inspection in post-op use, open design leaves most common access sites for Distal Radius Fracture surgery uncovered

## Adjustable

- Provides customized fit with high-comfort by allowing the HCP to adjust the aluminium core

## Versatile

- Ideal for managing fractures from the acute phase and stable, non-displaced fractures to post-op care, post-splinting and post-cast stabilization

# Benefits for your patients

## Supports Pain Relief

- 3-strap system with aluminum splint offers targeted stability while limiting wrist movements

## Hygienic

- Wrist brace is water-resistant and can be worn during hand washing and bathing

## Comfortable

- High-comfort and discreet wear due to the low-profile, open and light-weight design

## Easy to Wear

- Open-frame design for simple application even with the non-dominant hand

## Safe

- Can be worn even in contact with food thanks to food-grade material



# Multiple fields of application: when to use Manus Air

Conservative treatment of stable,  
non-displaced distal radius fractures

Post-operative stabilization of plated  
distal radius/ulna fractures

Post-traumatic stabilization  
after cast removal

Wrist sprains  
Tendinitis  
Carpal Tunnel Syndrome (CTS)  
Oeseoarthritis or rheumatoid arthritis  
(arthritis of the wrist)



## Place your order now

We'd love to talk to you about a trial order. Speak to your sales representative to find out more.

Actimove  
Catalog



Description	REF #	Measurement	Size
Actimove Manus Air Right, Small	6500500	6.8 cm – 8 cm / 2¾" – 3"	R, Small
Actimove Manus Air Right, Medium	6500501	8 cm – 9 cm / 3" – 3½"	R, Medium
Actimove Manus Air Right, Large	6500502	9 – 10 cm / 3½" – 4"	R, Large
Actimove Manus Air Left, Small	6500503	6.8 cm – 8 cm / 2¾" – 3"	L, Small
Actimove Manus Air Left, Medium	6500504	8 cm – 9 cm / 3" – 3½"	L, Medium
Actimove Manus Air Left, Large	6500505	9 – 10 cm / 3½" – 4"	L, Large

\*This product is not made with natural rubber latex.



Actimove®,  
an Essity brand



essity.com | actimove.com

BSN Medical Inc., an Essity company  
5825 Carnegie Blvd., Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325  
To order toll-free: Tel. (+1) 800 537 1063

65503 R1 ©2025 BSN Medical Inc. G25

The BSN Medical Inc., Reimbursement Disclaimer: The Centers for Medicare & Medicaid Services (CMS) uses Palmetto GBA LLC., as the Pricing, Data Analysis, and Coding (PDAC) Contractor. The PDAC Contractor maintains the Product Classification List, located at <https://www4.palmettogba.com/pdac.dmecls/>. The Product Classification List is a searchable database containing products that have received HCPCS coding verification from the PDAC and their corresponding HCPCS Codes. BSN Medical Inc., however, makes no representations as to the accuracy of the information contained within the PDAC Database, nor any representations as to whether its products are reimbursable under any government sponsored healthcare program and/or private-insurance program. It is the provider's sole responsibility, in consultation with the insurer, if necessary, to determine medical necessity, ensure coverage criteria is met, and submit appropriate codes, modifiers, and charges for services rendered. RV: 12/24