

# Actimove®

Feel good. Move better.

## *Rhizo Forte*

Stability and function  
for everyday activities




# Rhizo Forte

The patented, unique thumb orthosis for everyday wear

- Helps to alleviate pain arising from osteoarthritis of the thumb.<sup>1</sup>
- Supports pain relief by placing the thumb in a functional position and immobilizing the carpometacarpal (CMC) and the first metacarpophalangeal (MCP1) joints.
- Improved design for optimized fit and comfort around the thumb.
- For everyday wear, as wrist and fingers are free to move while the thumb is effectively restricted.
- Provides excellent fit and high comfort by customizing anatomical shape.
- Individually adjustable to the palm by molding the aluminum core.
- Maximum wearing comfort thanks to the open, low-profile design.
- Measure the palm width for the correct size.
- Easy to close with one hand using hook & loop straps.
- Easy cleansing of the brace by using hand wash.
- For easy daily hygiene: Water-resistant brace dries quickly after hand washing.
- This orthosis can be worn when in contact with food.<sup>2</sup>
- Suitable for those sensitive to latex, as this product is not made with natural rubber latex.
- Protected by patent, see [www.essity.com/patents](http://www.essity.com/patents).
- Available in Small, Medium and Large.
- Available in gray and black.
- For stock keeping, long shelf life of five years.

## Fields of application include:

- Thumb osteoarthritis
- Rheumatoid arthritis of the thumb
- Ulnar collateral ligament injuries of the MCP 1 joint (e.g., skier's thumb)
- Instability of the CMC joint
- Soft tissue injuries (sprains)
- Traumatic thumb dislocation
- Ligament instabilities
- Postsurgical and nonsurgical treatment

	Code		Size	Measurements	Color	Unit
	Left	Right				
	7623801	7623800	S	2¾" – 3" 6.8 – 8 cm	Gray	1 Thumb Orthosis
	7623803	7623802	M	3" – 3½" 8 – 9 cm	Gray	1 Thumb Orthosis
	7623805	7623804	L	3½" – 4" 9 – 10 cm	Gray	1 Thumb Orthosis
	7623819	7623818	S	2¾" – 3" 6.8 – 8 cm	Black	1 Thumb Orthosis
	7623821	7623820	M	3" – 3½" 8 – 9 cm	Black	1 Thumb Orthosis
	7623823	7623822	L	3½" – 4" 9 – 10 cm	Black	1 Thumb Orthosis



<sup>1</sup>Essity/BSN medical GmbH, 2021, data on file, submitted for publication.

<sup>2</sup>SGS group, 2021, data on file.



Actimove®,  
an Essity brand

 [essity.com](http://essity.com) | [actimove.com](http://actimove.com)

BSN Medical Inc., an Essity company  
5825 Carnegie Blvd., Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325  
To order toll-free: Tel. (+1) 800 537 1063