

# Push<sup>®</sup> med Back Brace

Supports the lower back and reduces pain


- Ensures proper positioning of the lumbar spine
- Individually adjusts each separate reinforcing element
- Adjusts compression via elastic bands, which increases the support of the metal stays
- Comfortable to wear during a wide range of activities
- Does not show through clothing
- Is machine washable on gentle cycle at low temperature
- Creates soft, comfortable feel and wicks moisture through the Sympress<sup>™</sup> microfibers

## Indications

- One vertebra slipping relative to another (spondylolisthesis)
- Crack or stress fracture of a vertebra (spondylolysis)
- Degenerative disc disease / osteoarthritis
- Lumbar compression fracture
- Post-operative support

## Measuring / Sizing

Measure around the waist as shown.



Please note that you should consult a trained professional for fitting and customization. BraceLab cannot accept brace returns that have been worn or altered.

SIZE	INCHES	CENTIMETERS
1	25 <sup>5</sup> / <sub>8</sub> - 29 <sup>1</sup> / <sub>2</sub>	65 - 75
2	29 <sup>1</sup> / <sub>2</sub> - 33 <sup>1</sup> / <sub>2</sub>	75 - 85
3	33 <sup>1</sup> / <sub>2</sub> - 38 <sup>1</sup> / <sub>8</sub>	85 - 97
4	38 <sup>1</sup> / <sub>8</sub> - 43 <sup>1</sup> / <sub>2</sub>	97 - 110
5	43 <sup>1</sup> / <sub>2</sub> - 49 <sup>1</sup> / <sub>4</sub>	110 - 125
6	49 <sup>1</sup> / <sub>4</sub> - 57 <sup>1</sup> / <sub>8</sub>	125 - 145

## Insurance Information

BraceLab is not a provider. We do not file insurance nor can we provide information on insurance coverage or reimbursement rates. Please contact your insurance provider for filing information.

