# Push® med Shoulder Brace Plus

## Shoulder immobilization without neck/shoulder strain

- Additional strap across the unaffected shoulder comfortably supports the hand and wrist, guaranteeing a stable position of the lower arm
- Adjustability of lower arm angle distributes arm weight comfortably across chest and shoulder straps, preventing blood flow constriction
- Is machine washable on gentle cycle at low temperature
- Creates soft, comfortable feel and wicks moisture through the Sympress™ microfibers

#### Indications

- Collar bone (clavicle) fracture
- Shoulder dislocation / subluxation
- Shoulder instability
- Post-operative shoulder care
- Frozen shoulder (adhesive capsulitis)
- Shoulder fracture (proximal humerus)

### Measuring / Sizing

Measure around the upper chest as shown.

Please note that you should consult a trained professional for fitting and customization. BraceLab cannot accept brace returns that have been worn or altered.

### **Insurance Information**

BraceLab is not a provider. We do not file insurance nor can we provide information on insurance coverage or reimbursement rates. Please contact your insurance provider for filing information.





#### Sizing: Right or Left Shoulder

SIZE	INCHES	CENTIMETERS
1	<b>21</b> <sup>5</sup> / <sub>8</sub> - <b>29</b> <sup>1</sup> / <sub>2</sub>	55 - 75
2	29 <sup>1</sup> / <sub>2</sub> - 41 <sup>3</sup> / <sub>8</sub>	75 - 105
3	41 <sup>3</sup> /8 - 55 <sup>1</sup> /8	105 - 140







