



ZERO-G Heel Cloud MODEL 9538, 9538W

The NYOrtho ZERO-G Heel Cloud™ combines craftsmanship and innovation providing a sides-down, open door design that provides an accessible, easy to monitor wound site. High-performance fabrics promote and maintain a cool and clean healing surface. Hypoallergenic microfiber filling comfortably cradles the foot, eliminating pressure and friction to the skin while offloading the heel and bony prominence on both sides of the ankle.



ANTIMICROBIAL



WATER
RESISTANT



COOL
BREATHABLE
FABRIC



DEVICE
FRIENDLY



OFF-LOADS

Wipe clean nylon exterior fabric is easy to maintain.

Soft interior peachskin fabric with texture gently helps maintain heel position minimizing patient movement.

Hypoallergenic microfiber filling comfortably cradles the foot and cushions the ankles, eliminating pressure while offloading the heel.

Breathable, water-resistant, antimicrobial peachskin fabric helps promote air circulation and wick moisture away.

Wedge stabilizer helps prevent foot and leg rotation.



Our patented **"fall away foot pad"** gives the practitioner the fullest and most unobstructed view of any heel suspension boot on the market.

THREE DVT HOLES

Left, right and center for attachment of DVT pressure devices.

SIDES-DOWN, OPEN DOOR DESIGN

Offers a direct, unobtrusive view of the wound without any contact or exerted pressure to the wound site.

OPEN BOTTOM OFFLOADS THE HEEL

Eliminates pressure and thereby pressure sores.

ELIMINATE PLANTAR FLEXION

Prevents plantar flexion with an easy to adjust positioner.