

CANDO cotton tensitube

- Create pressure garments or wrist and ankle braces
- Use liner under orthotic devices to help prevent hypertrophic scarring and contracture
- Applies pressure to limbs or trunk where injury is localized
- Vary compression by using either a single or double layer
- Cotton material
- Requires no pins or tape


Width	11 yards
2"	24-4040
2.5"	24-4041
2.75"	24-4042
3"	24-4043
3.5"	24-4044
4"	24-4045
4.75"	24-4046
6.75"	24-4047
8"	24-4048
10"	24-4049
12.75"	24-4050



CANDO stockinette

- Use liner under fracture casts and braces
- Conforms to the shape of the body
- Prevents wrinkles from forming under cast or brace
- **Cotton material** works great with plaster casts
- **Black polyester nylon material** adds durability to artificial limbs and other appliances
- Each roll is 25 yards in length

Width	Cotton	Black Polyester
1"	24-4053	24-4155
2"	24-4054	24-4156
3"	24-4055	24-4157
4"	24-4056	24-4158



black polyester stockinette

cotton stockinette

CANDO stockinette liners

- Precut liner pieces are fast and easy-to-use
- Use liner under hand and wrist splints, braces or casts
- Cotton material provides a smooth, wrinkle-free fit
- One size fits most patients
- Comes in a pack of 10
- Available in 2 colors, with or without thumb material

Color	with Thumbs	without Thumbs
natural	24-4160	24-4162
black	24-4161	24-4163



natural liner w/o thumbs

black liner w/ thumbs

CANDO heel / elbow protector

- Slip-on protector has layer of foam that cushions pressure points in either the knee or elbow
- Conforms to body and promotes air circulation
- Available in 4 sizes

Size	
small	24-4166
medium	24-4167
large	24-4168
x-large	24-4169



FABRICATION ENTERPRISES INC
 PO Box 1500, White Plains, NY 10602 (USA)
 tel: 1-800-431-2830 / 1-914-345-9300
 fax: 1-800-634-5370 / 1-914-345-9800
 FabEnt.com sales@FabEnt.com